

Rapid Labs

Diagnostic Product List 2016

www.rapidlabs.co.uk

Rapid Labs Ltd, Unit 2 & 2a Hall Farm Business Centre
Church Road, Little Bentley, Colchester
Essex. CO7 8SD, United Kingdom.

email: info@rapidlabs.co.uk or medical@rapidlabs.co.uk

Quality Diagnostic Products

Rapid Labs is an independent medical diagnostics business established in 2008, engaged in the production, sourcing and worldwide distribution of quality diagnostic reagents, test kits, dehydrated culture media, EIA, rapid tests, glass vials, glass bottles and dropper/pipette assemblies.

Our UK facilities are ISO 9001 & 13485 certified. We offer expertise that ensures that our raw material, production process, quality and costs are carefully controlled, production activities are streamlined and in accordance with international standards, whilst our finished products are quickly and efficiently dispatched and processed through customs with internationally accepted documentation.

The warehousing, manufacture, packaging, sales, distribution and administration functions of the business are conducted from our premises in the United Kingdom.

The Rapid Labs product range INDEX is classified into:

- ABO Blood Grouping reagents Page 3
- Lectin Blood Grouping reagents Page 4
- Rhesus Blood Grouping reagents Page 4
- Rare Blood Grouping reagents Page 4
- Blood Bank Support Products Page 4
- Rapid Latex Tests Page 5
- Syphilis Screening and Confirmation Page 5
- Febrile Antigens Page 6
- Elisa Kits Page 7 - 9
- Rapid Test Device Page 10 - 11
- Rapid Test Strip Page 12
- Drugs of Abuse Kits (Urine and Saliva Samples) Page 13
- Urinalysis Strips Page 14
- Bacteriology, Discs Rings, Supplements, Bacterial Identification. Page 14
- Chromogenic Dehydrated Culture Media Page 15
- Dehydrated Culture Media Page 15 – 18
- Notes Page 19
- Glass vials and bottles Page 20
- Plastic bottles and packaging Page 20
- *Available to order: chemicals, stains and dyes. Bacterial Typing Antisera for Meningitides, Salmonella, Vibrio Cholera and Shigella.*

Rapid Labs is a forward thinking organisation that responds quickly to customer demand creating unique solutions to our customer and market demands.



URS is a member of Registrar of Standards (Holdings) Ltd.




BLOOD GROUPING REAGENTS

ABOD, AHG, SEROLOGICAL ALBUMIN



Cat. No.	Product Description		Packing	Notes
BG-A5	Anti-A	Monoclonal	5 ml	NEW
BG-A10	Anti-A	Monoclonal	10 ml	
BG-B5	Anti-B	Monoclonal	5 ml	NEW
BG-B10	Anti-B	Monoclonal	10 ml	
BG-AB5	Anti AB (Clone 1)	Monoclonal	5 ml	NEW
BG-AB10	Anti AB (Clone 1)	Monoclonal	10 ml	
BG-AB5C2	Anti AB (Clone 2)	Monoclonal	5 ml	NEW
BG-AB10C2	Anti AB (Clone 2)	Monoclonal	10 ml	
BG-DBLEND5	Anti-D (IgG & IgM)	Monoclonal	5 ml	NEW
BG-DBLEND10	Anti-D (IgG & IgM)	Monoclonal	10 ml	
BG-DBHT10	Anti-D (IgG & IgM high titre)	Monoclonal	10 ml	NEW
BG-DIGM5	Anti-D IgM	Monoclonal	5 ml	NEW
BG-DIGM10	Anti-D IgM	Monoclonal	10 ml	
BG-DIGMHT10	Anti-D IgM high titre	Monoclonal	10 ml	NEW
BG-DIGMHT5	Anti-D IgM high titre	Monoclonal	5 ml	NEW
BG-DGG5	Anti D IgG	Monoclonal	5 ml	
BG-DGG10	Anti D IgG	Monoclonal	10 ml	
BG-ABD10	Anti A, B, D (IgG & IgM) kit	Monoclonal	3x10 ml	
BG-ABOD10	Anti A, B, AB, D (IgG & IgM) kit	Monoclonal	4x10 ml	
ANTI-HUMAN GLOBULIN & SEROLOGICAL ALBUMIN				
BG-AHG5	Anti-Human Globulin (Coombs) Green	Polyspecific	5 ml	NEW
BG-AHG10	Anti-Human Globulin (Coombs) Green	Polyspecific	10 ml	
BG-AHGM5	Anti-Human Globulin (Coombs) Green	Monospecific	5 ml	NEW
BG-AHGM10	Anti-Human Globulin (Coombs) Green	Monospecific	10 ml	
BG-AHGC5	Anti-Human Globulin (Coombs) Clear	Polyspecific	5ml	
BG-AHGC10	Anti-Human Globulin (Coombs) Clear	Polyspecific	10ml	
BG-BSA225	Bovine Albumin 22%		5 ml	NEW
BG-BSA2210	Bovine Albumin 22%		10 ml	
BG-BSA305	Bovine Albumin 30%		5 ml	NEW
BG-BSA3010	Bovine Albumin 30%		10 ml	

Bulk 1L pack sizes available for all Blood Grouping reagents





BLOOD GROUPING REAGENTS

RH-Hr SYSTEM REAGENTS,
RARE TYPING REAGENTS & LECTINS



Cat. No.	Product Description		Packing	Notes
BG-A1L5	Anti A1	Lectin	5 ml	
BG-A1L10	Anti-A1	Lectin	10 ml	
BG-H2	Anti H	Lectin	2 ml	
BG-H5	Anti-H	Lectin	5 ml	
BG-H10	Anti-H	Lectin	10 ml	
BG-NL2	Anti-N	Lectin	2 ml	
BG-NL5	Anti-N	Lectin	5 ml	
BG-BC5	Anti-C (big)	Monoclonal	5 ml	
BG-BC10	Anti-C (big)	Monoclonal	10 ml	
BG-LC5	Anti-c (little)	Monoclonal	5 ml	
BG-LC10	Anti-c (little)	Monoclonal	10 ml	
BG-CDE2	Anti-CDE	Monoclonal	2 ml	NEW
BG-CDE5	Anti-CDE	Monoclonal	5 ml	
BG-CDE10	Anti-CDE	Monoclonal	10 ml	
BG-BE5	Anti-E (big)	Monoclonal	5 ml	
BG-BE10	Anti-E (big)	Monoclonal	10 ml	
BG-LE5	Anti-e (little)	Monoclonal	5 ml	
BG-LE10	Anti-e (little)	Monoclonal	10 ml	
BG-KELL5	Anti-K (Kell)	Monoclonal	5 ml	
BG-KELL10	Anti-K (Kell)	Monoclonal	10 ml	
BG-CW2	Anti-Cw (Human)	Monoclonal	2 ml	
BG-RHCH10	Rh Control		10 ml	
POTENTIATORS				
BG-LISSWS10	Low Ionic Strength Solution (Working		10 ml	
BG-LISS10	Low Ionic Strength Solution + Add		10 ml	
RARE TYPING REAGENTS				
BG-P1M2	Anti-P1(Monoclonal)		2 ml	
BG-BS2	Anti-S (Monoclonal)		2 ml	
BG-SMSH2	Anti-s (Human)		2 ml	
BG-SMS2	Anti-s (Monoclonal)		2 ml	
BG-FYA2	Anti-Fya (Monoclonal)		2 ml	
BG-FYBH2	Anti-Fyb (Human)		2 ml	
BG-SMKH2	Anti-k (Human)		2 ml	
BG-KPAH2	Anti-Kpa (Human)		2 ml	
BG-KPBH2	Anti-Kpb (Human)		2 ml	
BG-JKAH2	Anti-Jka (Human)		2 ml	
BG-JKBH2	Anti -Jkb (Human)		2 ml	
BG-LUAH2	Anti -Lua (Human)		2 ml	
BG-LUBH2	Anti-Lub (Human)		2 ml	
BG-LEA2	Anti-Lea (Monoclonal)		2 ml	
BG-LEB2	Anti-Leb (Monoclonal)		2 ml	
BG-MH2	Anti M (Human)		2 ml	
BG-PAP10	Papenzyme (Papain)		10 ml	
BG-BROM10	Bromelite (Bromelin)		10 ml	
PR121	3 Cell Screen in Alsevers		3x10ml	NEW
PR144	Panel 1, 10 Cell Screen		10 x 5ml	NEW

Bulk 1L pack sizes available for all Blood Grouping reagents



LATEX AGGLUTINATION

50-100 TEST KITS AND BULK



Cat. No.	Product Description	Packing	Notes
RL-ASO100	ASO Latex Test Kit	100T/kit	
RL-ASO50	ASO Latex Test Kit	50T/kit	
RL-CRP100	CRP Latex Test Kit	100T/kit	
RL-CRP50	CRP Latex Test Kit	50T/kit	
RL-RA100	RA (RF) Latex Test Kit	100T/kit	
RL-RA50	RA (RF) Latex Test Kit	50T/kit	
RL-PREG100	Pregnancy Latex (indirect) Test Kit	100T/kit	
RL-PREG50	Pregnancy Latex (indirect) Test Kit	50T/kit	
RL-LE100	Systemic Lupus Erythematosus Latex Test Kit	100T/kit	
RL-LE50	Systemic Lupus Erythematosus Latex Test Kit	50T/kit	
RL-LE25	Systemic Lupus Erythematosus Latex Test Kit	25T/kit	
RL-MONO50	Mononucleosis Latex Test Kit	50T/kit	
RL-MONO20	Mononucleosis Latex Test Kit	20T/kit	
RL-STA50	Staphylococcus Latex Test Kit	50T/kit	
RL-RB100	Rose Bengal Agglutination Test Kit	100T/kit	NEW
RL-TOX100	Toxoplasmosis Latex Test Kit	100T/kit	NEW
RL-WR100	Waler Rose Haemagglutination	100T/kit	NEW

Bulk 1L pack size available for all Latex reagents



SYPHILIS TESTING

**PARTICLE AGGLUTINATION, TPHA, RPR & VDRL
Carbon Antigen (see also EIA & RAPIDS)**

Cat. No.	Product Description	Packing	Notes
RPR & VDRL			
D-RPR100	RPR Kit- 100T	100T/kit	
D-RPR250	RPR Kit- 250T	250T/kit	
D-RPR500	RPR Kit- 500T	500T/kit	
D-RPR500A	RPR Kit- 500T (no accessories)	500T/kit	NEW
D-RPRCARDS100	RPR Test Cards 10 well	10x10 cards	
RL-VDRL250	VDRL Carbon Antigen 5ml+ controls (no	5 ml + 2 x 1ml	NEW
RL-RPR5ML	RPR Carbon Antigen	5 ml	
RL-RPRP1ML	RPR Carbon Antigen Positive Control	1 ml	
RL-RPRN1ML	RPR Carbon Antigen Negative Control	1 ml	
TPHA PARTICLE AGGLUTINATION			
RL-TPHA100	TPHA	100T Kit	
RL-TPHA200	TPHA	200T Kit	
RL-TPHA500	TPHA	500T Kit	
RL-TPHA2000	TPHA for PK7200 and PK7300 analysers	2000T Kit	NEW
RL-SPC1	TPHA Positive Control	1ml	
RL-SNC1	Syphilitic Negative Control	1ml	

Bulk 1L pack size available for all Syphilis reagents



FEBRILE ANTIGENS

WIDAL KITS & CONTROLS



Cat. No.	Product Description	Packing	Notes
HIGH DENSITY FEBRILE ANTIGENS			
HF-BAA5	Brucella Abortus Antigen	5 ml	
HF-BMA5	Brucella Melitensis Antigen	5 ml	
HF-POXK5	Proteus OXK Antigen	5 ml	
HF-POX25	Proteus OX2 Antigen	5 ml	
HF-POX195	Proteus OX19 Antigen	5 ml	
HF-SOA5	Salmonella O Antigen Group A	5 ml	
HF-SOB5	Salmonella O Antigen Group B	5 ml	
HF-SOC5	Salmonella O Antigen Group C	5 ml	
HF-SOD5	Salmonella O Antigen Group D (Typhi O)	5 ml	
HF-SHA5	Salmonella H Antigen a	5 ml	
HF-SHB5	Salmonella H Antigen b	5 ml	
HF-SHC5	Salmonella H Antigen c	5 ml	
HF-SHD5	Salmonella H Antigen d (Typhi H)	5 ml	
F-FNC1	Febrile negative control sera	1 ml	
HF-FPC1	Febrile positive control sera	1 ml	
HF-KIT8X5	Widal Kit Salmonella Ha,Hb,Hc,Hd,OA,OB,OC,OD.	8 x 5ml	
HF-KIT8X5C	Widal Kit Salmonella Ha,Hb,Hc,Hd,OA,OB,OC,OD. +tve/-tve	8 x 5ml, 2 x 1ml	
STANDARD FEBRILE ANTIGENS			
F-BAA5	Brucella Abortus Antigen	5 ml	
F-BMA5	Brucella Melitensis Antigen	5 ml	
F-POXK5	Proteus OXK Antigen	5 ml	
F-POX25	Proteus OX2 Antigen	5 ml	
F-POX195	Proteus OX19 Antigen	5 ml	
F-SOA5	Salmonella O Antigen Group A	5 ml	
F-SOB5	Salmonella O Antigen Group B	5 ml	
F-SOC5	Salmonella O Antigen Group C	5 ml	
F-SOD5	Salmonella O Antigen Group D (Typhi O)	5 ml	
F-SHA5	Salmonella H Antigen a	5 ml	
F-SHB5	Salmonella H Antigen b	5 ml	
F-SHC5	Salmonella H Antigen c	5 ml	
F-SHD5	Salmonella H Antigen d (Typhi H)	5 ml	
F-FPC1	Febrile positive control sera	1 ml	
F-FNC1	Febrile negative control sera	1 ml	
F-KIT8X5	Widal Kit Salmonella Ha,Hb,Hc,Hd,OA,OB,OC,OD.	8 x 5ml	
F-KIT8X5C	Widal Kit Salmonella Ha,Hb,Hc,Hd,OA,OB,OC,OD. +tve/-tve	8 x 5ml + 2 x 1ml	

Bulk 1L pack size available for all Febrile Antigen reagents



ELISA & EIA TEST



Cat. No.	Product Description	Total incubation time (min)*	Sensitivity	Packing	Notes
SYPHILIS & CHAGAS					
DA-SYTA96	EIA- Syphilis Total Antibody	60 min	0.0015Iu/ml	96T	
DA-SYTA96	EIA- Syphilis Total Antibody	60 min	0.0015Iu/ml	96T	
DA-SYTA192	EIA- Syphilis Total Antibody	60 min	0.0015Iu/ml	96T	
DA-SYTA480	EIA- Syphilis Total Antibody	60 min	0.0015Iu/ml	192T	
DA-CHAG96	EIA- Chagas IgG	60 min	100%	480T	
FEVER					
DA-CLG96	Chlamydia IgG	30+30+20 min		96T	<i>NEW</i>
DA-CLM96	Chlamydia IgM	30+30+20 min		96T	<i>NEW</i>
DA-LEPM96	Leptospira IgM ELISA Kit			96T	
DA-DGG96	Dengue IgG ELISA Kit			96T	
DA-DAG96	Dengue Ag ELISA Kit			96T	
DA-DGM96	Dengue IgM ELISA Kit			96T	
DA-CHKGM96	CHIK IgM ELISA Kit			96T	
STEROID HORMONES (SERUM)					
SEK-0001	Cortisol	60+15 min	2,44 ng/mL	96T	
SEK-0002	Testosterone	60+15 min	0,07 ng/mL	96T	
SEK-0003	Estradiol	120+30 min	15 ng/mL	96T	
SEK-0004	17OH-Progesterone	60+15 min	0,1 ng/mL	96T	
SEK-0005	DHEA-S	60+15 min	0,045	96T	
SEK-0006	Progesterone	60+15 min	0,05 ng/mL	96T	
SEK-0007	Free Estriol	60+15 min	0,05 ng/mL	96T	
SEK-0008	Androstenedione	60+15 min	0,04 ng/mL	96T	
SEK-0015	Free Testosterone	60+15 min	0,10 pg/mL	96T	
SEK-0018	Urinary Cortisol	60+15 min	2 ng/mL	96T	
SEK-0019	Total Estriol	60+15 min	1 ng/mL	96T	
SEK-0053	Aldosterone	60+20 min	15 pg/mL	96T	
SEK-0068	Estrone			96T	
FERTILITY HORMONES					
SEK-0009	LH	60+15 min	0,22 IU/mL	96T	
SEK-0010	FSH	60+15 min	0,17 IU/mL	96T	
SEK-0011	Prolactin	60+15 min	0,12 ng/mL	96T	
SEK-0087	SHBG	60+15 min	0,2nmol/m	96T	<i>NEW</i>
SEK-10205	HCG	60+20 min	2.5 IU/l	96T	<i>NEW</i>
SEK-10235	Free beta HCG	30+30+20 min	2 ng/ml	96T	<i>NEW</i>
SEK-10238	PAPP-A	60+20 min	20 mu/l	96T	<i>NEW</i>
HEART AND LUNG DISEASES					
SEK-10250	CRP Ultra	30+30+20 min	0.05 mg/l	96T	<i>NEW</i>
SEK-10291	Troponin	60+20 min	0.25 ng/ml	96T	<i>NEW</i>
SEK-10238	Alveomucin	30+30+20 min	20 U/ml	96T	<i>NEW</i>
MISCALANIOUS					
SEK-0039	Ferritin	60+10 min	1 ng/mL	96T	
SEK-0050	hGH	60+15 min	0,8 µIU/mL	96T	

Cat. No.	Product Description	Total incubation time (min)*	Sensitivity	Packing	Notes
THYROID HORMONES					
SEK-10201N	TSH NEO	1080+60+20 min	3 µIU/ml	96T	<i>NEW</i>
SEK-0013	TSH	60+15 min	0,01 IU/mL	96T	
SEK-0037	FREE T3	60+15 min	0,05 pg/mL	96T	
SEK-0038	FREE T4	60+15 min	0,05 ng/dL	96T	
SEK-0044	T3	60+15 min	5 ng/dL	96T	
SEK-0045	T4	60+15 min	0,4 ug/dL	96T	
SEK-0048	THYROGLOBULIN	90+60+15 min	0,44 ng/mL	96T	
SEK-10131	aTPO	30+30+20 min	5 IU/ml	96T	<i>NEW</i>
SEK-10132	aTG			96T	<i>NEW</i>
TUMOR MARKERS					
SEK-0012	AFP	60+15 min	1 ng/ML	96T	
SEK-0051	CEA	60+15 min	1 ng/ML	96T	
SEK-0052	PAP	30+10 min	0,5 ng/mL	96T	
SEK-0054	CA125	60+15 min	1 IU/mL	96T	
SEK-0055	CA 15-3	60+60+15 min	5 U/mL	96T	
SEK-0056	CA 19-9	60+15+15 min	5 U/mL	96T	
SEK-0057	PSA Total	60+15 min	0,5ng/mL	96T	
SEK-0064	FREE PSA	60+15 min	0,052ng/m	96T	
SEK-0073	h-NSE new	60+15 min	0,075	96T	
SEK-0074	S100B	120+60+20 min	10 pg/mL	96T	
SEK-10244	CA 72-4	120+60+20 min	2 U/ml	96T	<i>NEW</i>
SEK-10227	MUC1 M22	30+30+20 min	1.5 U/ml	96T	<i>NEW</i>
SEK-10228	MUC1 M20	60+20 min	1 U/ml	96T	<i>NEW</i>
DIABETES					
SEK-0076	Insulin	120+15 min	2 µIU/mL	96T	
SEK-0077	C-Peptide new	120+15 min	0,025	96T	
SEK-0082	GAD	60+15+15 min	92,3%	96T	
SEK-0083	IIA (Anti Insulin)	60+60+20+20 min	60%	96T	
SEK-0084	IA2 (Anti (protein/tyrosine	O.N.+60+20+20 min	75%	96T	
AUTOIMMUNE					
SEK-0095	anti-DsDNA IgG	30+30+15 min	93.30%	96T	
SEK-0098	anti-ENA Screen	30+30+15 min	94,7%	96T	
SEK-0099	anti-ANA Screen	30+30+15 min	83,9%	96T	
SEK-10133	aMPO	30+30+20 min	3 IU/ml	96T	<i>NEW</i>
SEK-10200	Total IgE	30+30+20 min	6 IU/ml	96T	<i>NEW</i>
INFECTIOUS					
SEK-10106	Mycoplasma IgG	30+30+20 min		96T	<i>NEW</i>
SEK-10121	Aspergillus IgG	30+30+20 min		96T	<i>NEW</i>
SEK-10126	Ureaplasma IgG			96T	<i>NEW</i>
SEK-10171	Giardia lamblia Total Ab			96T	<i>NEW</i>

Minimum order quantity 15 Pieces of each kit



COMPACT ELISA TEST KITS



Cat. No.	Product Description	Total incubation time (min)*	Sensitivity	Packing	Notes
INFECTIOUS DISEASES					
E-101196	HIV-1/2 O	60 min	>99.9%	96T	Qual.
E-128196	HIV-1/2 O Antigen/Antibody	120 min	>99.9%	96T	Qual.
E-129196	HTLV I/II	90 min	>99.9%	96T	Qual.
E-102196	HBsAg	70 min (60/10)	0.2ng/ml*	96T	Qual.
E-103196	HCV	70 min	>99.9%	96T	Qual.
E-105196	HBsAb	45 min (30/15)	99.3%	96T	Qual.
E-127196	HbsAb Quantitative	45 min (30/15)	99.3%	96T	Quant.
E-106196	HBeAg	45 min (30/15)	99.7%	96T	Qual.
E-107196	HBeAb	45 min (30/15)	99.1%	96T	Qual.
E-108196	HBCAb TA (IgG+IgM+IgA)	45 min (30/15)	99.1%	96T	Qual.
E-121196	HEV IgM	75 min	99.2%	96T	Qual.
E-123196	H pylori Antigen (Ag) (incl.extraction kit)	70 min (60/10)		96T	Qual./Quant.
E-124196	H pylori IgG	70 min		96T	Qual./Quant.
E-131196	HBCAb IgG	70min	99.1%	96T	Qual.
TORCH					
E-109196	Toxo IgG	70 min	>99.9%	96T	Qual./Quant.
E-110196	Toxo IgM	70 min	>99.9%	96T	Qual.
E-111196	Rubella IgG	70 min	96.4%	96T	Qual./Quant.
E-112196	Rubella IgM	70 min	90.0%	96T	Qual.
E-113196	CMV IgG	70 min	>99.9%	96T	Qual./Quant.
E-114196	CMV IgM	70 min	89.3%	96T	Qual.
E-115196	HSV 1 IgG	70 min	98.5%	96T	Qual.
E-116196	HSV1 IgM	70 min	96.3%	96T	Qual
E-117196	HSV 2 IgG	70 min	92.1%	96T	Qual.
E-118196	HSV 2 IgM	70 min	90.5%	96T	Qual.
E-119196	HSV 1/2 IgG	70 min	98.6%	96T	Qual.
E-120196	HSV 1/2 IgM	70 min	85.7%	96T	Qual.
THYROID					
E-301196	TSH	75 min (60/15)	0.078uIU/mL	96T	Quant.
E-304196	Total T3	75 min (60/15)	0.1 ng/mL	96T	Quant.
E-302196	Total T4	75 min (60/15)	0.4 ug/dL	96T	Quant.
E-305196	Free T3	75 min (60/15)	0.24 pg/mL	96T	Quant.
E-303196	Free T4	75 min (60/15)	0.06 ng/dL	96T	Quant.
FERTILITY					
E-401196	Rapid-hCG	15 min (10/5)	0.6 mIU/mL	96T	semi-Quant.
E-402196	LH	75 min (60/15)	0.5 mIU/mL	96T	Quant.
E-404196	Prolactin	75 min (60/15)	1.0 ng/mL	96T	Quant.
E-403196	FSH	75 min (60/15)	0.5 mIU/mL	96T	Quant.
TUMOR					
E-201196	AFP (Alpha-Fetoprotein)	45 min (30/15)	0.6 ng/mL	96T	Quant.
E-202196	CEA	45 min (30/15)	0.38 ng/mL	96T	Quant.
E-203196	PSA Total	45 min (30/15)	0.1 ng/mL	96T	Quant.



RAPID TEST DEVICES

10 - 40 TEST INDIVIDUALLY POUCHED DEVICES



Cat. No.	Product Description	Size	Notes
FERTILITY			
D-HCGD20	HCG Test - Urine/serum/plasma	Devices 20T	
D-HCGD40	HCG Test - Urine/serum/plasma	Devices 40T	
D-LHD20	Ovulation (LH) Test - Urine	Devices 20T	
RAPID - INFECTIOUS DISEASE TEST PRODUCTS			
D-HAVMD20	HAV IgM Rapid Test	Devices 20T	
D-HBSGD20	HBsAg Test - whole blood/serum/plasma	Devices 20T	
D-HBSGD40	HBsAg Test - whole blood/serum/plasma	Devices 40T	
D-HBSGSPD40	HBsAg Test - serum/plasma	Devices 40T	NEW
D-HBSBD20	HBsAb Test- whole blood/serum/plasma	Devices 20T	
D-HBEGD20	HBeAg Test -whole blood/serum/plasma	Devices 20T	
D-HBEBD20	HBeAb Test - whole blood/serum/plasma	Devices 20T	
D-HBCBD20	HBcAb Test - whole blood/serum/plasma	Devices 20T	
D-HBCMD20	HBcIgM Test - whole blood/serum/plasma	Devices 20T	
D-HEVMD20	HEV IgM Rapid Test	Devices 20T	
D-HBVM5D10	Multi HBV-5 Test - whole blood/serum/plasma	Multi device 10T	
D-HCVD20	HCV Test - whole blood/serum/plasma	Devices 20T	
D-HCVD40	HCV Test - whole blood/serum/plasma	Devices 40T	
D-HCVSPD40	HCV TEST – serum/plasma	Devices 40T	NEW
D-HIV0WBD20	HIV 1,2,O tri-line Test - whole blood/serum/plasma	Devices 20T	
D-HIV0WBD40	HIV 1,2,O tri-line Test - whole blood/serum/plasma	Devices 40T	
D-HIV0SPD40	HIV 1,2,O tri-line Test - serum/plasma	Devices 40T	NEW
D-HIV1-2D20	HIV 1/2 Ab Plus Combo - whole blood/serum/plasma	Devices 20T	
D-HIVAGABD20	HIV Ag/Ab 4th Gen. - whole blood/serum/plasma	Devices 20T	NEW
D-HIVAGABD40	HIV Ag/Ab 4th Gen. - whole blood/serum/plasma	Devices 40T	NEW
D-HIVSYPD20	HIV/Syphilis Ab -serum/plasma	Devices 20T	
D-TPD20	TP(syphilis) Test - whole blood/serum/plasma	Devices 20T	
D-TPD40	TP(syphilis) Test - whole blood/serum/plasma	Devices 40T	
D-TPSPD40	TP(syphilis) Test - serum/plasma	Devices 40T	NEW
D-SCREENMD10	Multi Infectious(HBsAg/HCV/HIV/TP)Test - whole blood/serum/plasma	Multi device 10T	
D-CHABD20	Chagas Ab Rapid Test whole blood/serum/plasma	Devices 20T	
SEXUALLY TRANSMITTED DISEASE RAPID DEVICE TEST			
D-CHLAGD20	Chlamydia Ag Sample SWAB	Devices 20T	
D-GOND20	Gonorrhea Sample SWAB	Devices 20T	
D-HPVD20	HPV Test	SWAB 20T	NEW
AUTOIMMUNE			
D-RFD20	RF Test whole blood/serum/plasma	Devices 20T	NEW
FEVER RAPID TESTS			
D-DGMD20	Dengue IgG/IgM Combo Test - whole blood/serum/plasma	Devices 20T	NEW
D-DAGMD20	Duo Dengue Ag-IgG/IgM	Devices 20T	NEW
D-DAGD20	Dengue Ag	Devices 20T	NEW
D-DGMCMD20	Dengue + CHIK (IgG/IgM-CHIK IgM)	Devices 20T	NEW
D-FILGMD20	Filariasis IgG/IgM Combo	Devices 20T	NEW
D-CHIKMD20	Chikungunya IgM whole blood/serum/plasma	Devices 20T	NEW
D-LEPGMD20	Leptospira IgG/IgM Combo	Devices 20T	NEW

Cat. No.	Product Description	Size	Notes
FEVER RAPID TESTS (CONTINUED)			
D-TBD20	TB Test - whole blood/serum/plasma	Devices 20T	
D-TOXMD20	TOXO IgM Test - whole blood/serum/plasma	Devices 20T	
D-TOXGD20	TOXO IgG Test - whole blood/serum/plasma	Devices 20T	
D-MONOD25	Mononucleosis Test - whole Blood/serum/plasma	Devices 25T	NEW
TORCH			
D-CMVMD10	CMV IgM Test- serum/plasma	Devices 10T	NEW
D-RUBMD10	Rubella IgM Test- serum/plasma	Devices 10T	NEW
D-HSV12D10	HSV ½ IgM Test- serum/plasma	Devices 10T	NEW
D-TORCHD25	ToRCH IgM Test- serum/plasma	Devices 25T	NEW
D-HSV1D20	HSV-1 IgG/IgM Test whole blood/serum/plasma	Devices 20T	NEW
D-HSV2D20	HSV-2 IgG/IgM Test whole blood/serum/plasma	Devices 20T	NEW
D-TOXOD10	TOXO IgG/IgM - serum/plasma	Devices 10T	NEW
PARASITIC RAPID TESTS			
D-MPFD20	Malaria Pf Ag - Whole Blood Device	Devices 20T	
D-MPFPVD20	Malaria Pf/Pv Ag Device	Devices 20T	
D-MPFPAND20	Malaria Pf/Pan Ag Device	Devices 20T	
D-MPFPVBD20	Malaria Pf/Pv Ab Combo	Devices 20T	
D-LEIGMD20	Leishmania IgG/IgM Combo	Devices 20T	
GASTRO INTESTINAL RAPID TESTS			
D-HPD20	H.Pylori Ab - whole blood/serum/plasma	Devices 20T	
D-HPAGD20	H. pylori Ag	Devices 20T	
D-CHAG20D	Cholera Ag	Devices 20T	
D-TYGM20D	Typhoid IgG/IgM serum/plasma	Devices 20T	
D-TYGMCD20	Typhoid IgG/IgM whole blood/serum/plasma	Devices 20T	
D-ROTAGD20	Rotavirus Ag Rapid Test	Devices 20T	
D-ROAAGD20	Rota/Adeno Ag Rapid Test	Devices 20T	
D-NOROD25	Norovirus Test- faeces	Devices 25T	NEW
RAPID - CANCER MARKER TEST PRODUCTS			
D-AFPD20	AFP Test - whole blood/serum/plasma	Devices 20T	
D-FOBD10	Faecal Occult Blood(FOB) Test - faeces	Devices 10T	
D-PSAWBD20	PSA Test - whole blood/serum/plasma	Devices 20T	
D-PSASPD20	PSA Test - serum/plasma	Devices 20T	
D-PSA2NGD20	PSA Test 2ng	Devices 20T	NEW
RAPID - CARDIAC MARKER TEST PRODUCTS			
D-TROPD20	Troponin I Test - whole Blood/serum/plasma	Devices 20T	
D-MCKTMD20	Multi-3(Myoglobin/CK-MB/Troponin I)panel - whole	Devices 20T	NEW
D-ACED40	Acetaminophen (ACE) Test - whole Blood/serum/plasma	Devices 40T	NEW
D-CTNLD10	cTnl Test - whole Blood/serum/plasma	Devices 10T	NEW
D-CRPD10	C-reactive protein (CRP) Test - whole Blood/serum/plasma	Devices 10T	NEW
D-DIMERD10	D-dimer Test - whole Blood/serum/plasma	Devices 10T	NEW
D-PCTD40	Procalcitonin (PCT) Test - serum/plasma	Devices 40T	NEW



RAPID TEST STRIPS

10-100 TEST INDIVIDUALLY POUCHED STRIPS



Cat. No.	Product Description	Size	Notes
RAPID - FERTILITY TEST PRODUCTS			
D-HCGUS25	HCG Test – Urine only	Canister Strips 25T	<i>NEW</i>
D-HCGUS50	HCG Test – Urine only	Strips 50T	<i>NEW</i>
D-HCGUS100	HCG Test – Urine only	Strips 100T	<i>NEW</i>
D-HCGS50	HCG Test - Urine/serum/plasma	Strips 50T	
D-HCGS100	HCG Test - Urine/serum/plasma	Strips 100T	
D-LHS50	Ovulation (LH) Test - Urine	Strips 50T	
RAPID - INFECTIOUS DISEASE TEST PRODUCTS			
D-HBSGS50	HBsAg Test - whole blood/serum/plasma	Strips 50T	
D-HBSGS100	HBsAg Test - whole blood/serum/plasma	Strips 100T	
D-HBSGSP50	HBsAg Test - serum/plasma	Strips 50T	<i>NEW</i>
D-HBSBS50	HBsAb Test- whole blood/serum/plasma	Strips 50T	
D-HBEGS50	HBeAg Test -whole blood/serum/plasma	Strips 50T	
D-HBEB50	HBeAb Test - whole blood/serum/plasma	Strips 50T	
D-HBCBS50	HBcAb Test - whole blood/serum/plasma	Strips 50T	
D-HBCMS50	HBcIgM Test - whole blood/serum/plasma	Strips 50T	
D-HCVS50	HCV Test - whole blood/serum/plasma	Strips 50T	
D-HCVS100	HCV Test - whole blood/serum/plasma	Strips 100T	
D-HCVSPS50	HCV Test – serum/plasma	Strips 50T	<i>NEW</i>
D-HIV0WBS50	HIV 1,2,O tri-line Test - whole blood/serum/plasma	Strips 50T	
D-HIV0WBS100	HIV 1,2,O tri-line Test - whole blood/serum/plasma	Strips 100T	
D-HIV0SPS50	HIV 1,2,O tri-line Test - serum/plasma	Strips 50T	
D-TPS50	TP (syphilis) Test - whole blood/serum/plasma	Strips 50T	
D-TPS100	TP (syphilis) Test - whole blood/serum/plasma	Strips 100T	
D-TPSPS50	TP (syphilis) Test - serum/plasma	Strips 50T	<i>NEW</i>
FEVER RAPID TEST			
D-TOXMS50	TOXO IgM Test - whole blood/serum/plasma	Strips 50T	
D-TOXGS50	TOXO IgG Test - whole blood/serum/plasma	Strips 50T	
D-DGMS50	Dengue IgG/IgM Tri Line WB/S/P Rapid Test	Strips 50T	
D-STEPS20	Strep A	SWAB 20T	
D-STRB20	Strep B	SWAB 25T	<i>NEW</i>
D-INFAS20	Influenza A/B	SWAB 20T	
RAPID – CANCER MARKER TEST PRODUCTS			
D-AFPS50	AFP Test - whole blood/serum/plasma	Strips 50T	
D-FOBS10	Faecal Occult Blood(FOB) Test - faeces	Strips 10T	
D-PSAWBS50	PSA Test - whole blood/serum/plasma	Strips 50T	
D-PSASPS50	PSA Test - serum/plasma	Strips 50T	
SEXUALLY TRANSMITTED DISEASE RAPID DEVICE TEST			
D-GONS50	Gonorrhea Rapid Test SWAB	Strips 50T	
PARASITIC RAPID TESTS			
D-LEIABS50	Leishmania Ab	Strips 50T	
GASTRO INTESTINAL RAPID TESTS			
D-HP50	H.Pylori Ab - whole blood/serum/plasma	Strips 50T	



RAPID DRUGS OF ABUSE

STRIPS AND DEVICES



Cat. No.	Product Description	Size	Notes
RAPID – DRUGS OF ABUSE TEST			
D-DOA1S50	Amphetamine (AMP) Test - Urine	Strips 50T	
D-DOA2S50	Methamphetamine (MET) Test -	Strips 50T	
D-DOA3S50	Opiates/Morphine (OPI) Test - Urine	Strips 50T	
D-DOA4S50	Barbiturates (BAR) Test - Urine	Strips 50T	
D-DOA5S50	Benzodiazepines (BZO) Test - Urine	Strips 50T	
D-DOA6S50	Cocaine (COC) Test - Urine	Strips 50T	
D-DOA7S50	Methadone (MTD) Test - Urine	Strips 50T	
D-DOA8S50	Cannabinoids (THC) Test - Urine	Strips 50T	
D-DOA9S50	Ketamine (KAN) Test - Urine	Strips 50T	
D-DOA10S50	Tricyclic Antidepressants (TCA) Test -	Strips 50T	
D-DOA11S50	Buprenorphine (BUP) Test - Urine	Strips 50T	
D-DOA12S50	Ecstasy (MDMA) Test - Urine	Strips 50T	
D-DOA13S50	Phencyclidine (PCP) Test - Urine	Strips 50T	
RAPID – DRUGS OF ABUSE DEVICES			
D-DOA1D20	Amphetamine (AMP) Test – Urine	Devices 20T	
D-DOA2D20	Methamphetamine (MET) Test -	Devices 20T	
D-DOA3D20	Opiates/Morphine (OPI) Test - Urine	Devices 20T	
D-DOA4D20	Barbiturates (BAR) Test - Urine	Devices 20T	
D-DOA5D20	Benzodiazepines (BZO) Test - Urine	Devices 20T	
D-DOA6D20	Cocaine (COC) Test - Urine	Devices 20T	
D-DOA7D20	Methadone (MTD) Test - Urine	Devices 20T	
D-DOA8D20	Cannabinoids (THC) Test - Urine	Devices 20T	
D-DOA9D20	Ketamine (KAN) Test - Urine	Devices 20T	
D-DOA10D20	Tricyclic Antidepressants (TCA) Test -	Devices 20T	
D-DOA11D20	Buprenorphine (BUP) Test - Urine	Devices 20T	
D-DOA12D20	Ecstasy (MDMA) Test - Urine	Devices 20T	
D-DOA13D20	Phencyclidine (PCP) Test - Urine	Devices 20T	
D-DOA14D20	Amphetamine (AMP) Test - Saliva	Devices 20T	
D-DOA15D20	Benzodiazepines (BZO) Test - Saliva	Devices 20T	
D-DOA16D20	Cannabinoids (THC) Test - Saliva	Devices 20T	
D-DOA17D20	Cocaine (COC) Test - Saliva	Devices 20T	
D-DOA18D20	Methamphetamine (MET) Test -	Devices 20T	
D-DOA19D20	Opiates/Morphine (OPI) Test - Saliva	Devices 20T	
D-DOA20D20	Ecstasy (MDMA) Test Cassette -	Devices 20T	
RAPID – DRUGS OF ABUSE MULTI-			
D-DOAM5U	DOA Multi 5 Parameter- Urine Test	Multi Device 10T	
D-DOAM10U	DOA Multi 10 Parameter- Urine Test	Multi Device 10T	
D-DOAM12U	DOA Multi 12 Parameter- Urine Test	Multi Device 10T	
D-DOACM12	DOA Multi 12 Parameter –Urine Cup	Multi Device 25T	NEW

***Please contact us for the latest list of available tests**



URINALYSIS STRIPS

TUBES OF 50-100 TEST STRIPS



Cat. No.	Product Description	Size	Notes
D-U2100	Urinalysis Strips 2 Parameter (Glucose, Protein)*	100T	
D-U3100	Urinalysis Strips 3 Parameter	100T	
D-U4100	Urinalysis Strips 4 Parameter	100T	
D-U5100	Urinalysis Strips 5 Parameter	100T	
D-U6100	Urinalysis Strips 6 Parameter	100T	
D-U7100	Urinalysis Strips 7 Parameter	100T	
D-U8100	Urinalysis Strips 8 Parameter	100T	
D-U9100	Urinalysis Strips 9 Parameter	100T	
D-U10100	Urinalysis Strips 10 Parameter	100T	
D-U11100	Urinalysis Strips 11 Parameter	100T	
	* Other combinations available		
D-U3410100	Urinalysis 10P Europe 100T Strips	100T	
D-U2MC100	Urinalysis Strips (Micro Albumin and Creatinine)	100T	
D-U2MC50	Urinalysis Strips (Micro Albumin and Creatinine)	50T	



ANTIBIOTIC DISCS & RINGS

IDENTIFICATION REAGENTS

Cat. No.	Product Description	Size	Notes
	ANTIBIOTIC SENSITIVITY DISCS: VIALS AND CARTRIDGES		
	All commercially utilised antibiotics available	5x 50 Discs Cartridges	
		Vial 50 or 100 discs	
	LYOPHILISED ANTIBIOTICS IN VIALS		
	All commercially utilised antibiotics available	Pack of 10	
	METHICILLIN RESISTANCE		
	Strip	50 Strips	
	MICROBIAL IDENTIFICATION DISCS		
	All Currently Utilised antibiotics available.	Sealed Tins of 6,4 or rings	
		3-Tipped 1x50	
	MICROBIAL IDENTIFICATION RINGS		
CM-CAT	Catalase Reagent	Dropper Bottle	
DD-BET-A	Betalactamase Strips	50 Strips	
DD-COAG-C	Coagulase Citrated Plasma	5 x 2ml vials	
DD-COAG-E	Coagulase EDTA Plasma	5 x 2ml vials	
	ISOLATION SUPPLEMENTS		
	All Currently Utilised antibiotics available.	Vials	Pack of 10



DEHYDRATED CULTURE MEDIA

500gm WRAPPED POTS



Cat. No.	Product Description	Size	Notes
CHROMOGENIC DEHYDRATED CULTURE MEDIA			
CM-CM01	Chr. Candida	500gm	
CM-CM02	Sakazaki GBP	500gm	
CM-CM03	E Coli	500gm	
CM-CM04	Listeria Aloa	500gm	
CM-CM05	Salmonella	500gm	
CM-CM06	Uricolour	500gm	
CM-CM07	Staphilocolor	500gm	
DEHYDRATED CULTURE MEDIA			
CM-AAB190	Aeromonas Agar Base	500gm	
CM-AB317	Agar Powder	500gm	
CM-AA318	Agar base	500gm	
CM-ATMC048	Amies Transport Medium (modified)w.charcoal (1%)	500gm	
CM-AIA151	Anaerobic Isolation Agar	500gm	
CM-APT370	APT agar	500gm	
CM-ADB299	Azide dextrose broth	500gm	
CM-BCM289	Bacillus Cereus Medium	500gm	
CM-BPA203	Baird Parker Agar	500gm	
CM-BHIA175	BHI Agar	500gm	
CM-BHIB106	BHI Broth	500gm	
CM-BEAA259	Bile Esculin azide agar	500gm	
CM-BSA283	Bismuth Sulphite Agar	500gm	
CM-BAB109	Blood Agar Base	500gm	
CM-BAB120	Blood Agar Base No.2.	500gm	
CM-BGA765	Bordet gengou agar	500gm	
CM-BTB053	Brewers Thioglycollate Broth	500gm	
CM-BGA182	Brilliant Green Agar	500gm	
CM-BGA167	Brilliant Green Agar Modified	500gm	
CM-BGBB127	Brilliant Green Bile Broth 2%	500gm	
CM-BAB131	Brucella Agar Base	500gm	
CM-BPW049	Buffered Peptone Water	500gm	
CM-CAB254	Campylobacter agar base	500gm	
CM-CBM130	Campylobacter Blood free Medium	500gm	
CM-CEB072	Campylobacter Enrichment Broth	500gm	
CM-CBTM054	Cary Blair transport Medium	500gm	
CM-CA148	Cetrimide Agar	500gm	
CM-CBLA261	China blue lactose agar	500gm	
CM-CLED096	CLED DI Medium	500gm	
CM-CLED097	CLED SI Medium	500gm	
CM-CAB137	Columbia Agar Base	500gm	
CM-CMM214	Cooked Meat medium	500gm	
CM-CMA159	Corn meal agar	500gm	
CM-DCAH173	DCA Agar Hynes	500gm	
CM-DCAL162	DCA Agar Leifsons	500gm	
CM-DTM0198	Dermatophyte Test Medium	500gm	
CM-DTA077	Dextrose Tryptone agar	500gm	
CM-DSSA708	Diagnostic Semi-Solid Salmonella Agar	500gm	
CM-DSA216	Dextrose Starch Agar	500gm	
CM-DNAA134	Dnase Agar	500gm	
CM-DST123	DSTA Agar	500gm	

Cat. No.	Product Description	Size	Notes
CM-DRCM156	DRCM medium	500gm	
CM-ECB263	EC broth	500gm	
CM-ECB250	EC broth modified	500gm	
CM-EEB161	EE Broth (Mossel)	500gm	
CM-EB355	Enterococcus broth	500gm	
CM-EA356	Enterococcus agar	500gm	
CM-EYA376	Egg Yolk Agar	500gm	
CM-EMB110	EMB agar	500gm	
CM-EA132	Endo Agar	500gm	
CM-FBB189	Fraser Broth Base	500gm	
CM-GCAB116	GC Agar Base	500gm	
CM-GCB401	Giolitti Cantoni broth	500gm	
CM-GPB372	Glucose Phosphate Broth	500gm	
CM-GNB387	Gram-negative broth	500gm	
CM-HEA207	Hektoen Enteric Agar	500gm	
CM-HMB210	Hoyles Medium Base	500gm	
CM-ISA031	Iron Sulphite Agar	500gm	
CM-ISA085	Iso-Sensitest agar	500gm	
CM-KEA154	Kanamycin-Esculin Agar	500gm	
CM-KEB089	Kanamycin-Esculin Broth	500gm	
CM-KIA194	Kligler Iron Agar	500gm	
CM-KCM195	Koser citrate medium	500gm	
CM-LB035	Lactose Broth	500gm	
CM-LPB258	Lactose phosphate broth	500gm	
CM-LBB266	LB Broth	500gm	
CM-LBA265	LB Agar	500gm	
CM-LTSB297	Lauryl Tryptose (Sulphate) Broth APHA	500gm	
CM-LB139	Lethen Broth	500gm	
CM-LEB098	Listeria Enrichment Broth FDA	500gm	
CM-LIMB188	Listeria Isolation Med. Base Oxford Form	500gm	
CM-LM236	Litmus Milk	500gm	
CM-LTMIND	Litmus powder indicator	100gm	
CM-LJM351	Löwenstein-Jensen Medium	500gm	
CM-LDB267	Lysine decarboxylase broth	500gm	
CM-LIA090	Lysine Iron Agar	500gm	
CM-MAWS185	MacConkey Agar (base)	500gm	
CM-MA2174	MacConkey agar No.2	500gm	
CM-MA3172	MacConkey Agar No.3	500gm	
CM-MAS170	MacConkey Agar Sorbitol	500gm	
CM-MAN257	MacConkey Agar w/o Crystal Violet	500gm	
CM-MANS158	MacConkey Agar w/o Salt	500gm	
CM-MBP094	MacConkey Broth Purple	500gm	
CM-MA381	Malt Agar	500gm	
CM-MEA164	Malt Extract Agar	500gm	
CM-MEB050	Malt Extract Broth	500gm	
CM-MALTEX	Malt Extract	500gm	
CM-MEGA169	Malt East 50% Glucose Agar	500gm	
CM-MSA213	Mannitol Salt Agar	500gm	
CM-MSB244	Mannitol Salt Broth	500gm	
CM-MRD029	Maximal recovery diluent	500gm	
CM-EMF387	m-Endo MF broth	500gm	
CM-MLS095	Membrane lauryl sulphate broth	500gm	
CM-MIO217	MIO medium	500gm	
CM-MIM111	Motility-Indol-Urea Base (MIU)	500gm	
CM-MB380	Motility Broth	500gm	
CM-MGM128	Modified Glutamate Medium	500gm	
CM-MRSA205	MRS Agar	500gm	
CM-MRSB177	MRS Broth	500gm	

Cat. No.	Product Description	Size	Notes
CM-MRVP272	MRVP medium	500gm	
CM-MSRV086	M.S.R.V	500gm	
CM-MSM041	MS Medium	500gm	
CM-MHA135	Mueller Hinton Agar NCCLS	500gm	
CM-MKB211	Mueller Kaufmann Broth	500gm	
CM-NB402	Nitrate broth	500gm	
CM-NA073	Nutrient Agar	500gm	
CM-NBE036	Nutrient Broth Economy	500gm	
CM-NG102	Nutrient Gelatine	500gm	
CM-NB2037	Nutrient Broth No.2 APHA	500gm	
CM-OGYE093	O.G.Y.E. Agar Base	500gm	
CM-PW038	Peptone Water	500gm	
CM-PB334B	Peptone bacteriological from bovine meat	500gm	
CM-PB252	Peptone broth	500gm	
CM-PAB150	Perfringens Agar Base	500gm	
CM-PAA032	Phenylalanine Agar	500gm	
CM-PRA108	Phenol Red Agar	500gm	
CM-PRB058	Phenol Red Broth	500gm	
CM-PRLB220	Phenol Red Lactose Broth	500gm	
CM-PRSB219	Phenol Red Sucrose Broth	500gm	
CM-PCAAP064	Plate Count Agar APHA	500gm	
CM-PCASP051	Plate Count Agar SP	500gm	
CM-PDA117	Potato Dextrose Agar	500gm	
CM-PDB118	Potato Dextrose Broth	500gm	
CM-PAB160	Pseudomonas Agar Base	500gm	
CM-PAB561	Purple Agar Base	500gm	
CM-R2AM043	R 2 A Medium	500gm	
CM-RVB070	Rapport Vassiliadis Broth	500gm	
CM-RCA168	Reinforced Clostridium Agar	500gm	
CM-RBAB083	Rose Bengal Agar Base	500gm	
CM-SDAC393	Sab Dex agar Chloramp 0,05 g/l	500gm	
CM-SDA197	Sab Dex Agar	500gm	
CM-SDB323	Sabouraud dextrose broth	500gm	
CM-SDAC200	Sab Dex agar Chloramphenicol 0,5 g/l	500gm	
CM-SLM157	Sabouraud Liquid Medium USP	500gm	
CM-SMA201	Sabouraud Maltose Agar	500gm	
CM-SSB244	Salmonella Shigella Broth	500gm	
CM-SA140	Schaedler Agar	500gm	
CM-SDAB255	SDA Beer Medium	500gm	
CM-SBB046	Selenite Broth Base	500gm	
CM-SCB047	Selenite-Cystine Broth	500gm	
CM-SMA375	Semisolid Media Agar	500gm	
CM-STA121	Sensitivity Test Agar	500gm	
CM-SIM081	SIM Medium	500gm	
CM-SCA063	Simmons Citrate Agar	500gm	
CM-SBA143	Slanetz Bartley Agar	500gm	
CM-SSA202	SS Agar	500gm	
CM-SSB249	SS Broth	500gm	
CM-SA055	Starch agar	500gm	
CM-STM042	Stuart Transport Medium	500gm	
CM-TCBS183	TCBS	500gm	
CM-TA384	Teepol Agar	500gm	
CM-TBB153	Tetrathionate Broth Base	500gm	
CM-TUFM069	Thioglycolat USP Fluid medium	500gm	
CM-THB104	Todd Hewitt Broth	500gm	
CM-TW039	Tryptone Water	500gm	
CM-T338	Tryptone	500gm	

Cat. No.	Product Description	Size	Notes
CM-TGYA221	Tryptone glucose Yeast Extract Agar	500gm	
CM-TPB078	Tryptose Phosphate Broth	500gm	
CM-TSA107	TSA	500gm	
CM-TSB080	TSB	500gm	
CM-TSCA369	TSC Agar	500gm	
CM-TSI204	TSI	500gm	
CM-UAB066	Urea Agar Base	500gm	
CM-UBB028	Urea Broth Base	500gm	
CM-VRBA114	Violet Red Bile Agar (lactose)	500gm	
CM-VRBG115	Violet Red Bile Glucose	500gm	
CM-VJA196	Vogel & Johnson Agar	500gm	
CM-WCA274	Wilkins-Chalgren agar	500gm	
CM-WLNA209	WL Nutrient Agar	500gm	
CM-WA178	Wort Agar	500gm	
CM-XLD179	XLD Agar	500gm	
CM-YA060	Yeast Agar	500gm	
CM-YEA059	Yeast Extract Agar	500gm	
CM-YSA192	Yersinia Selective Agar (CIN)	500gm	
CM-YNB119	Yeast Nitrogen Base Broth	500gm	
CM-YEXPOWD	Yeast Powder	500gm	

Notes:

REQUEST OUR GLASS CATALOGUE FOR:

Glass dropper vials

Pipette/Dropper assemblies

Glass moulded bottles

HDPE bottles

PET bottles and vials

Plain caps and closures

Tamper evident caps

Atomisers

and much more... please see our website for details:

www.rapidlabs.co.uk