

ENDOSCOPY STORAGE CABINETS

Felmed ENT Scope Storage Cabinets offer effective drying and safe storage of flexible endoscopes. They are certified to the latest EN Standards, using the latest technologies for tracking and traceability, and offer 99.995% efficient HEPA Filtration

Certified to BS EN/16442 2015

Meets latest specifications

Ultra quiet fan with rapid air purge when doors are opened

Easy to clean, all stainless steel construction

User-friendly touch-sensitive control panel

Clear interlocked doors

Bar code reader and printer for recordable

tracking and traceability

Heavy duty foot master castors for stability and mobility

Ergonomic, easy to use scope supports, with

dedicated individual easy to clean compartments

Minimal cross contamination

Ultra-clean air to 0.15 microns

Optimised efficient drying

Available in 3 models; 5 scopes, 10 Scopes and 20 Scopes



FELCON

Rapid Labs



Height adjustable and static stainless steel sinks

Felmed manufacture a range of height adjustable and static stainless steel sinks – HTM64 compliant single or multiple height adjustable sink units for manual cleaning of both instruments and scopes. Optional washing spray guns can be supplied. Bespoke stainless steel height adjustable sinks are designed and manufactured in our factory to suit specific applications.

We work closely with our customers to ensure that we produce quality manufactured height adjustable sinks to fit exactly into the available space.

Stainless steel height adjustable electrically operated table

Robust, sturdy design in 304 or 316 stainless steel, with push button operation, 4 programmable user heights incorporating 'True Parallel Drive', which ensures that the tabletop will always remain level and true during use.

Tables can be supplied in any width, depth, height and load carrying range to suit the customer's individual requirements, with mountings for fan/filter modules, power points, equipment or storage racks which can be fitted to the worktop as required.

Authorised Distributor - **Rapid Labs Ltd**, Unit 1, 2 & 2a Hall Farm Business Centre, Church Road
Little Bentley, Colchester, Essex, CO7 8SD, UK

sales@rapidlabs.co.uk. +44 (0)1206 250484 www.rapidlabs.co.uk